

THE  FRANKENMAN ENDOLC PORTFOLIO
ENDOSCOPIC LINEAR CUTTING STAPLERS



 **FRANKENMAN**

AN UNPARALLELED BALANCE OF QUALITY AND COST EFFECTIVENESS



WHEN YOU PUT PATIENTS' HEALTH FIRST,
REMARKABLE THINGS HAPPEN




- 1 LARGER JAW OPENING
- 2 RELOADABLE CARTRIDGES
- 3 3 AVAILABLE SHAFT LENGTHS
- 4 ARTICULATION CONTROL
- 5 RETURN KNOB
- 6 ROTATION CONTROL
- 7 SAFETY SWITCH
- 8 GRANULATED SCALE

FRANKENMAN

OTHERS




COMPARISON OF JAW OPENING




AVAILABLE IN STRAIGHT OR ARTICULATED CARTRIDGES WITH UP TO 45 DEGREE BI-DIRECTIONAL ARTICULATION


ARTICULATING RELOADS

 30MM ARTICULATING LOADING UNIT
AVAILABLE IN WHITE, BLUE & GREEN



 45MM ARTICULATING LOADING UNIT
AVAILABLE IN WHITE, BLUE & GREEN



 60MM ARTICULATING LOADING UNIT
AVAILABLE IN WHITE, BLUE & GREEN



STRAIGHT RELOADS



30MM STRAIGHT LOADING UNIT

AVAILABLE IN WHITE, BLUE & GREEN



45MM STRAIGHT LOADING UNIT

AVAILABLE IN WHITE, BLUE & GREEN



60MM STRAIGHT LOADING UNIT

AVAILABLE IN WHITE, BLUE & GREEN



STAPLER DEVICE	CODE	SHAFT LENGTH (mm)	WORKING LENGTH (mm)	EXTERNAL DIAMETER (mm)	PER BOX
ENDOLC SHORT	ELC50	50	280	12.5	3
ENDOLC STANDARD	ELC150	150	380		
ENDOLC XL	ELC250	250	480		



ORDERING CODES AND SPECIFICATIONS

	CODE	OPEN STAPLE HEIGHT (mm)	CLOSED STAPLE HEIGHT (mm)	STAPLE ROWS	STAPLE LINE LENGTH (mm)	STAPLE QUANTITY	STAPLE CROWN	PER BOX
STRAIGHT	○ ESLCR 30W	2.5	1.0	6	33	48	3.0MM	6
	○ ESLCR 45W	2.5	1.0	6	45	66		
	● ESLCR 45B	3.8	1.5	6				
	● ESLCR 45Y	4.2	1.75	6				
	● ESLCR 45G	4.5	2.0	6				
	○ ESLCR 60W	2.5	1.0	6	61	90		
	● ESLCR 60B	3.8	1.5	6				
	● ESLCR 60Y	4.2	1.75	6				
	● ESLCR 60G	4.5	2.0	6				
ARTICULATING	○ ECLCR 30W	2.5	1.0	6	33	48		
	● ECLCR 30B	3.8	1.5	6	45	66		
	○ ECLCR 45W	2.5	1.0	6				
	● ECLCR 45B	3.8	1.5	6				
	● ECLCR 45Y	4.2	1.75	6				
	● ECLCR 45G	4.5	2.0	6	61	90		
	○ ECLCR 60W	2.5	1.0	6				
	● ECLCR 60B	3.8	1.5	6				
	● ECLCR 60Y	4.2	1.75	6				
	● ECLCR 60G	4.5	2.0	6				

ALSO AVAILABLE FROM FRANKENMAN



Frankenman are proud to offer several training courses via the Frankenman Academy, in numerous surgical disciplines, with world renowned expert faculty.

These courses are staged to cover both didactic elements and "hands on" skills workshops.

QUALITY CONTROL

We pride ourselves on world leading continual quality improvement processes.

ALL stapling devices and reloads have a unique serial number, as well a batch/lot number, to allow for more accurate tracking and traceability.

In addition we test fire 100% of our products (before they are then reloaded and packaged), to verify correct B-Shaped closure.

