

THE  FRANKENMAN
SUCTION & IRRIGATION SYSTEM





 FRANKENMAN




WHEN YOU PUT PATIENTS' HEALTH FIRST,
REMARKABLE THINGS HAPPEN

THE FRANKENMAN SUCTION AND IRRIGATION SYSTEM






-  The Suction and Irrigation system from FRANKENMAN offers unrivalled quality and unrivalled specifications in Suction and Irrigation... when you put patients' health first remarkable products follow.
-  The FRANKENMAN design creates a system focused on performance. High Capacity tubing coupled to friction reducing double "o" ring sealed pistons for both Irrigation and Aspiration have created a system with unique fluid flow characteristics... maximising flow rates... virtually eliminating clogging.



ALTERNATIVE GRIP CONFIGURATION



KEY FEATURES

-  MAXIMUM CONTROL - ERGONOMICALLY DESIGNED PISTOL GRIP
-  MAXIMUM VERSATILITY - CONVERTS TO ALTERNATIVE "PALM", INLINE GRIP CONFIGURATION
-  MAXIMUM APPLICATIONS - 3MM, 5MM & 0MM CANNULAS IN 280MM, 320MM & 420MM WORKING LENGTHS
-  INTEGRAL LUER LOCK TISSUE TRAP CONNECTOR
-  MAXIMUM QUALITY - MANUFACTURED TO THE HIGHEST PROFESSIONAL STANDARDS



ORDERING CODES AND SPECIFICATIONS

PRODUCT CODE	CONFIGURATION	CANNULA		PER BOX
		LENGTH	DIAMETER	
SIS328	TUBING SUPPLIED - WITH DOUBLE SPIKE CONNECTIONS	280MM	3MM	10
SIS528			5MM	
SIS1028			10MM	
SIS332		320MM	3MM	
SIS532			5MM	
SIS1032			10MM	
SIS342		420MM	3MM	
SIS542			5MM	
SIS1042			10MM	
SIS328T	DEVICE AND CANNULA ONLY	280MM	3MM	
SIS528T			5MM	
SIS1028T			10MM	
SIS342T		420MM	3MM	
SIS542T			5MM	
SIS1042T			10MM	

ALSO AVAILABLE FROM FRANKENMAN

STANDARD CIRCULAR STAPLERS



COMPACT CIRCULAR STAPLERS



LINEAR STAPLERS



LINEAR CUTTING STAPLERS



TRANSCUT™ - TRANSVERSE LINEAR CUTTING STAPLERS




CIRCULAR STAPLERS FOR RECTAL PROLAPSE & ANOPEXY



UNIVERSAL TRANSANAL ACCESS PLATFORM



ENDOSCOPIC LINEAR CUTTER



ACCESSORIES



frankenman.com/academy

Frankenman are proud to offer several training courses via the Frankenman Academy, in numerous surgical disciplines, with world renowned expert faculty.


These courses are staged to cover both didactic elements and "hands on" skills workshops.

QUALITY CONTROL

We pride ourselves on world leading continual quality improvement processes.

ALL stapling devices and reloads have a unique serial number, as well a batch/lot number, to allow for more accurate tracking and traceability.

In addition we test fire 100% of our products (before they are then reloaded and packaged), to verify correct B-Shaped closure.



www.frankenman.com | info@frankenman.co.uk | +86-512-68786150

Frankenman International Ltd - 108 Jinfeng South Road, High New District, - Suzhou, 215163, P.R. China